



Adelphi's Apollo range of Liquid Viewers offer a simple yet effective method of detecting particulate matter present in liquids.

Apollo I Liquid Inspection Unit

Designed for use in pharmaceutical and laboratory environments, the Apollo I Liquid Viewer has a polarised inspection light for the detection of particulate matter in transparent liquids contained in ampoules, vials and bottles. Particulate contamination can be unintentionally present in injections and intravenous infusions, and this unit provides a simple method of detection.

To operate, simply hold the container between the light panel and the polarising screen, shake or swirl the container, and any particles will show as bright specks against a dark background.



Key Features:

- Compact and hygienic stainless steel body
- Adjustable viewing angle to suit individual working preferences
- Adjustable magnifier
- LED light source for optimum bulb life
- Tool free assembly
- Offers a clear, polarised and magnified image
- Includes padded storage bag for safe storage & transport




Item code	Specification
6031301	APOLLO I LIQUID VIEWER
Length (mm)	320
Width (mm)	295
Height (mm)	420 (max)
Weight (kg)	Nett: 4.9 Gross: 7.2
Carton size (cm):	50 x 41 x 27
Electricity:	110/240V AC 50/60HZ

Apollo II Liquid Inspection Unit (for transparent glass)

Designed for use in pharmaceutical and laboratory environments, the Apollo II Liquid Viewer provides a simple method for the detection of visible particles in ampoules, vials, and bottles.

The unit is operated by simply swirling each individual container in front of both panels, and observing for any particulate contamination (other than air bubbles) unintentionally present in the solution. In front of the black panel, light coloured particulates will be visible, and in front of the white panel, dark coloured particulates will be visible.

This unit conforms to the European Pharmacopoeia Specification 2.9.20 and to the US Pharmacopoeia – USP chapter 41, section 6.1.



Key Features:

- Adjustable lamp unit for comfortable working height
- Uses two energy efficient LED lamps
- Fold up frame for space saving and easy storage - the inspection surfaces fold together to reduce the chance of damage
- Removable diffuser for use with coloured glass
- Conforms to European Pharmacopoeia 2.9.20 and to the US Pharmacopoeia – USP chapter 41, section 6.1
- Intuitive assembly with no tools required

Item code	Specification
6031511	APOLLO II LIQUID VIEWER
Length (mm)	714
Width (mm)	420
Height (mm)	450
Weight (kg)	Nett: 10 Gross: 7.9
Carton size (cm):	76 x 52.5 x 26.5
Electricity:	100/240V AC 50/60HZ



Supporting IQ/OQ documentation can be supplied upon request.

User manuals are available in English, Czech, Dutch, French, German and Spanish – additional languages are also available upon request.

Apollo II High Intensity Liquid Inspection Unit (for coloured glass)

The Apollo II High Intensity specification is ideal for coloured glass containers. It uses an updated three LED bulbs, giving a light level of 8,000 – 10,000 lux.

This unit also conforms to the European Pharmacopoeia Specification 2.9.20 and to the US Pharmacopoeia – USP chapter 41, section 6.1.



Supporting IQ/OQ documentation can be supplied upon request.





**Adelphi
Manufacturing**

A personal approach, on a global scale

www.adelphi.uk.com